

1.5 inches

Table 3.7 Parameters of Laser Used in Patient no.3 (Lesion No.2) During the Three Treatment Sessions

Spacing 24 lines

1 inch

Parameter	Session								
	1	2	3	4	5	6	7	8	9
Duration (μ s)	2000	2000	1600	1200	1200	1500	3000	3100	3200
Duration (mJ/cm ²)	23.8	23.8	19	14.3	14.3	17.9	35.7	36.9	38.1
Repeat (s)	2	1	single	1	1	single	1	1	single
Overlap (TH)	1	1	1	2	1	1	1	1	1
Distance (mm)	0.8	1	1	1	1	1	1	1	1
Shape	square	square	circle	square	square	circle	square	square	circle
Dots	420	420	216	420	420	216	420	420	217

1.5 inches

1 inch

Page No.

1 inch